




OXITEST

Oxidation Stability Reactor

Oxidative Stability Studies Directly on the Whole Sample
According to AOCS Standard Procedure



OXITEST Oxidation Stability Reactor

The OXITEST is designed for R&D and Quality Control laboratories interested in investigating the behaviour of fat oxidation stability. The OXITEST speeds up the oxidation process because of the two accelerating factors: temperature and pressure. The instrument measures the absolute pressure change of oxygen inside the two titanium heating chambers at constant temperature, thus monitoring the oxygen uptake by reactive components in the sample and automatically generating an IP value for further analysis via software.

REPRESENTATIVE RESULTS

The OXITEST ensures representative results by determining the oxidation stability of the whole sample without extensive sample preparation: be it solid, liquid, or doughy.

PREMIUM RESISTANCE

The oxidation chambers, sample holders and covers are made of titanium, a premium material guaranteeing:

- High resistance and excellent chemical compatibility;
- Cost savings by eliminating required consumables.

EASE OF USE

The OXITEST is entirely controlled via PC by the powerful OXISoft™ Software.

- All information at a glance;
- Multilingual support;
- Preloaded methods for data analysis.


PRODUCTIVITY

The OXITEST has 2 independent oxidation chambers and further expansion by connecting up to 4 systems to the same PC allowing the simultaneous analysis of up to 8 samples.



AOCS Standard Procedure Cd 12c-16

Accelerated Oxidation Test for the Determination of the Oxidation Stability of Foods, Oils and Fats using the Oxitest Oxidation Test Reactor



INDUCTION PERIOD - IP

Results

Induction Period (IP)	14h 46 min (Graphical method)
Test duration	21h 00 min
Curve 1	$Y = -0,003X + 6,18$
Curve 2	$Y = -1,575X + 29,43$

The test allows the sample oxidation curve to be obtained, characterized by an **Induction Period (IP)**. The Induction Period is the time required to reach the starting point of oxidation, corresponding to either a level of detectable rancidity or a sudden change in the rate of oxidation. **The longer the Induction Period, the higher the stability against oxidation over time.** The operator can create **test reports** for a single test or compare different analyses for a better interpretation of the data.

OXISoft™ Software

The OXISoft™ is available in different languages and comes with a pre-installed library of methods related to a wide range of sample types. The operator can use and modify them, or create personalized methods. In order to obtain a visible oxidation flex, the sample tested should contain 2-4% of unsaturated fatty acids. When product degradation is mainly due to the loss of the aromatic components and the oxidation flex is not visible, information on the product oxidation can be achieved by interpreting the slope of the oxidation curve. Many parameters can be investigated, including:




1 REPEATABILITY TEST

A series of tests run on the same sample or standard to verify its IP, to calculate accuracy and repeatability of the data.




2 FRESHNESS TEST

To verify the quality of different lots, for example of the same raw material, by comparing them. This can be valuable to confirm, for example, whether the product freshness is related to the cost of raw materials.




3 FORMULAS COMPARISON

Which ingredients are required to create the most stable formula of a finished product, under the same conditions. By comparing the results, OXISoft™ will be able to automatically distinguish the best formula, easily recognizable by the higher IP.




4 PACKAGING COMPARISON

Particularly useful for testing which packaging maintains the product in the freshest condition.



5 IP DURING AGEING


This procedure allows to obtain a prediction of the oxidation stability during a shelf life study by measuring the product at defined time intervals.



6 ESTIMATED SHELF LIFE TEST

It is possible to have a prediction of oxidation stability for shelf life studies. By following a dedicated procedure and testing the same product at different temperatures, in the case of a linear equation, the operator can extrapolate and estimate the oxidation stability of the sample even at room temperature.

VELP ERMES CONNECTION




Connect the OXITEST to the exclusive VELP Ermes Cloud Platform to improve your laboratory experience. The VELP Ermes Cloud Platform connection will unburden you from tedious tasks, improving your lab productivity and enhancing service support.


- Real-time monitoring and control of the instrument from PC, smartphone and tablet whenever you want, wherever you are;
- Immediate alert and notification with the possibility to stop the instrument for maximum safety;
- Regular software updates guaranteeing the best performance and new features with just one-click.

ermes enabled

RESISTANT TITANIUM PARTS



OXISoft™ POWERFUL AND INTUITIVE SOFTWARE



ANALYSIS ON THE WHOLE SAMPLE

INSTRUMENT - CODE

OXITEST	230 V / 50-60 Hz	F30900248
OXITEST	115 V / 60 Hz	F30910248

SUPPLIED WITH



10002948
OXITEST OXISoft™
Software



10003134
USB cable




10001985*
Sample holder




10001984*
Spacer



A00000236
High temperature
sealing greaser



10002033
O-ring 3475
viton



E00010012
VELP Ermes
1 Year Connection

* OXITEST is supplied with 6 sample holders and 4 spacers.



OPTIONAL ACCESSORIES

OXITEST IQ/OQ Manual	A00000242
Calibration temperature kit for Oxitest	A00000360
Calibration temperature kit for Oxitest without thermometer and probe	A00000376
High temperature sealing grease	A00000236
VELP Ermes 1 Year Connection	E00010012
VELP Ermes 3 Years Connection	E00010036

FIELDS OF APPLICATION

OXITEST works directly on the whole sample without the need for preliminary fat separation, ensuring representative results on solid, semi-solid and liquid samples, raw and finished products



FOOD, FEED AND BEVERAGE INDUSTRY



COSMETIC INDUSTRY



PHARMACEUTICAL INDUSTRY

TECHNICAL DATA

	OXTITEST
NUMBER OF OXIDATION CHAMBERS	2
CAPACITY OF SINGLE CHAMBER	up to 100 ml
MAX DEVIATION FROM THE SET TEMPERATURE	≤ 0.5 °C
REPRODUCIBILITY OF SET TEMPERATURE	≤ ± 0.2 °C
INTERFACE	USB
CONNECTIVITY	Cloud via LAN or Wi-Fi
POWER	900 W
POWER SUPPLY	115 V / 60 Hz - 230 V / 50-60 Hz
WEIGHT	16.5 Kg 36.3 lb
DIMENSION (WxHxD)	365x190x485 mm 14.6x7.6x19.4 in
OVERPRESSURE	Safety valve
OUT-RANGE TEMPERATURE	Visual alarm
DAMAGED PROBE	Visual alarm
TEMPERATURE RANGE	From room temp. to 120 °C
PRESSURE RANGE	0 - 8 bar



SERVICE & SUPPORT

VELP Scientifica products are designed by our engineers to resist years of laboratory use.

Our products are manufactured with premium materials to guarantee the best performance with maximum safety.

According to our experience, a proper and regular maintenance is necessary to ensure the highest performance of analytical instrument. VELP Service Department and VELP Official Partners are always ready to offer you maintenance and service support tailored to your needs.

GET THE SUPPORT YOU NEED CHOOSING THE OPTIONS:

- Installation
- Preventive Maintenance
- Help-desk and Remote support
- Technical Assistance
- Analytical Support
- Calibration Certification



We reserve the right to make technical alterations
We do not assume liability for errors in printing, typing or transmission

VELP Official Partner

LAT GmbH
Labor- und Analysen-Technik
Heinkelstr. 4
DE - 30827 Garbsen



DESIGNED AND MANUFACTURED IN ITALY



ITALY - HQ
Via Stazione 16
20865 Usmate (MB) Italy
Tel. +39 039 628811
velpitalia@velp.com

INDIA
velpindia@velp.com

USA
155 Keyland Court, Bohemia
NY 11716 - U.S.
Tel. +1 631 573 6002
velpusa@velp.com

CHINA
Xinlong Rd Building 28, Lane 1333
Shanghai city - China
Tel. +86 180 1755 7329
velpchina@velp.com



ermes

CONNECT
YOUR PRODUCT
NOW

SER 158 Series

www.velp.com

How to Connect to VELP Ermes Platform

SER 158 Series

In order to speed up the connection, be sure to in contact with your IT specialist to get all your laboratory intranet connection details.

Requirements:

1. **Device:** PC, Tablet, Smartphone with your intranet (LAN or Wi-Fi)
2. **Browser:** Chrome, Safari, Firefox, Edge, UC browser

FIRST STEP | ACTIVATE VELP ERMES ACCOUNT


A

Go to www.velp.com through your internet browser.



B

In the Smart Lab menu click on “Configure Your VELP Ermes Account”. Fill in the form and accept Terms and Conditions.



C

You will receive an email containing your ID, Password and a link to proceed to VELP Ermes authentication at the email address that you entered during the registration.



D

At the first login on VELP Ermes platform a new password will be requested: remember to change it!



SECOND STEP | CONNECT TO YOUR SER 158

Connect the LAN cable to the SER 158 using the LAN connection on the unit.

A

In the ControlPad Select Set-Up → Connectivity and set the IP address, subnetmask and gateway of your intranet.

N.B: Contact your IT specialist to get all your laboratory intranet connection details.

B

In the ControlPad Select Set-Up → Connectivity and press the down arrow until you reach **ermes** and select "product registration".

C

Select **YES** under the box "Enable SER 158 to transmit on the ERMES platform". **Do not exit from this page, keep it open.**

D

Login on VELP Ermes with your ID and Password received by email, and select "Device Registration".

E

Enter the SER 158 serial number, product code (you can find these information on the silver label on the back of it) and purchase date. Press Next.

F

Check the ControlPad where you will find the One Time Code (OTC) of four digits. Enter this number into VELP Ermes platform. Press Next and wait for connection.

G

Give a name to your SER 158 and enter the address where your instrument is located. Press Next.

H

Assign your SER 158 to a laboratory or create a new one selecting "Laboratory Management".

I

Now from "Instruments" you can monitor and control your connected SER 158.



Always **Connected** and **Easy to Use**

Monitor and manage multiple instruments 24/7 from your device. Real time visualization of your analysis and of the instrument working conditions. Immediate notifications and alerts.



Safely Manage **Your Data**

VELP Ermes collects and stores your data with the maximum level of encryption following the highest cybersecurity standards. You are now able to have access to your data anytime, anywhere. VELP Ermes will allow you to easily generate and share reports of your analyses.



Your Instrument **Always Updated**

Get software updates (OTA) immediately and benefit of all new features coming from VELP continuous product innovation.



Enhanced **Service Support**

Faster and smarter Service Support thanks to the direct contact with VELP specialists. With VELP Ermes we can guide you to get the maximum from your instrument.