

The last generation microbiological safety laminar flow cabinets **CRUMA BIO CLASS II Type A2** with digital functions.

The cabinet design ensures complete integrity by close control of the airflow. Incoming air is filtered through the first HEPA filter and down through the work area in a laminar flow pattern. The descending air creates a protective barrier preventing any outside air entering the cabinet. The airflow is then re-directed from the base of the cabinet into a plenum where 30% is expelled through a second HEPA filter and 70%, along with 30% new make up air, is re-circulated back into the cabinet through the first HEPA filter.



APPLICATIONS

- √ Manipulation of microorganisms, bacteria, fungi, viruses and parasites Risk categories 1, 2 and 3.
- √ Isolation and sample culture
- √ Quantification methods
- \checkmark Microscopy techniques and sample preparation Identification and classification of microorganisms
- √ Genetic Manipulation

MAIN SPECIFICATIONS

- $\checkmark\,$ Microprocessor controlled motor blower (LCD) with volumetric sensor for exhausted air flow monitoring
- \checkmark State of the art Microprocessor control system offering:
 - Large screen monitor
 - Automatic control of preset airflow volumes
 - Sliding sash window with smart control
 - Permanent monitoring of HEPA filters life span
 - Alarms. Multilevel alarms, with redundancy functions Permanent display of working conditions
 - Highest air flow stability both in case of transitional disturbances or to progressive filter clogging
 - Semi-automatic fumigation cycle (EN12297 tested and certified)

More information www.cruma.es

- Continuous monitoring of front barrier air flow for the highest operator safety
- Low barrier alarm
- Power failure alarm
- \checkmark Volt-free contact for remote monitoring of exhaust fan
- \checkmark Automatic reset of initial conditions in case of powerfailure
- \checkmark C-shaped support stand for the easiest one man installation procedure

CHARACTERISTICS

- √ Fully EN12469 certified by TÜV Hamburg
- \checkmark State of the art microprocessor control system
- \checkmark Large digital display, high resolution
- ✓ Air and aerosol-tight sliding sash, electrically operated by finger touch
- $\checkmark~$ Alarms for low airflow and wrong front window position
- \checkmark Sloped front and back wall for the most comfortable access
- \checkmark Front access for filter maintenance and service
- ✓ C-shaped support stand for the easiest one man installation procedure
- √ Easy retrofit option kits

TECHNICAL SPECIFICATIONS					
Number of filters HEPA-H14		2			
Number of IP44 fans		1			
Average volume of treated air		±350 m³/h			
Average face velocity		0.50 m/s			
Total electrical power consumption		260 W			
Voltage-Frequency		110-230 V / 50-60 Hz			
Fluorescent Lamp / Light intensity		36 W / 900 Lux			
Noise level		49 dB			
Packaging: Palletized cardboard box	Volume	2 m ³			
	Weight	260 Kg			

Do you need help or *technical assistance*?

Contact your distributor or call us if you have any questions or need technical support, spare parts,

maintenance service...

(R) our engagement

+34 93 370 61 62

2 year *warranty*

Because we are convinced of the quality of our products.





We recognise our responsability and dependence towards a healthy environment and, therefore, we destinate more than 7% of our annual budget in innovating and developing new products for the lab operator

cour 🔁 engagement

cruma.es export@cruma.es

CRUMABIO-1



SIZES (mm)					
External Dimensions		Internal Dimensions			
Widht 1075	Depth 840	Height 1450	Widht Depth Heig 924 600 70		

MAIN STRUCTURE	
External metal parts	1.2 mm galvanized steel, coated with antiacid polymer powder resin thermo-hardened at 200 °C
Internal metal parts	Stainless Steel internal surfaces with 2B finishing (including spillage tray). Solid work surface in 3 sections and special designed front grill at the air barrier that cannot be obstructed by the operator's arms
Front door	Electrically operated from control panel sliding multilayer safety glass window

ACCORDING TO STANDARDS			
Cabinet	CEN-12469		
Filters	EN-779: 1996 (HEPA & ULPA Filters) EN-1822:1998 (HEPA & ULPA Filters)		
Quality	UNE EN ISO 9001:2015		



KEYBOARD MEMBRANE

Membrane keyboard commands for opening and closing the front window, activation of service outlets, electrical socket, gas solenoid valve, lighting and predisposed key to activate the connector kit germicidal lamp under the condition of complete closure of front of the machine, given the presence of the relevant position sensors.

All available accessories can be mounted without any structural modification of the cabin.

Provision for possible expulsion for channel adapters or additional filters.







MECHANICAL AND FUNCTIONAL SPECIFICATIONS

Sloped front design for the highest operational comfort. Sloped back side of the working chamber for the best down flow disbution (cabinet carcass EN12298 tested and certified for air tightness)

Utilities inlets from the top of the cabinet

Stainless Steel internal surfaces with 2B finishing (including spillage tray). Solid work surface (3 sections) and special designed front grill

Electrically operated sliding multilayer safety glass window (max opening at 120°)

Comfortable 200mm front opening

Easy to install retrofit options through lateral sides.

Exposed exhaust HEPA filter for easy visual integrity check

H14 class High Efficiency Particulate Airfilters with 99.995% efficiency on.3micron particles (most penetrating particle diameter) (EN 1822-1 and EN 13091:1999 tested and certified)

Both exhaust and Main Filters are equipped with a micromesh membran elocated downstream which acts as airspeed equalizer expansion plenum, as well as a clear indicator of filter damages

Filter change and maintenance from the front of the cabinet

Exhaust transitions easily installable

Key operated. The key can be removed when the unit is in SAFE mode, in order to avoid unwanted operation. In case of power failure, the cabinet is reset to original working conditions

Self calibration cycle performed when cabinet is switched on

High speed rinse and set up cycle performed, before reaching the SAFE operating mode.

Visual display of SAFE conditions. Pre-warning before actual alarm conditions are reached (visual and acoustic alarms)

Soft touch control with keys for standard service utilities. Interconnected UV and fluorescent lights

Exhaust and recirculating flow rates ensure 25 air changes/min in the working area (30% 70% split)

Front barrier air speed ≥0.5mt/sec

Aperture protection Factor (Apf)≥1.5 x 10 exp5

Cleanability Index CC grade. (EN 12296 tested and certified)

SERIAL EQUIPMENT

DOP test port 2 Electrical sockets

UV lamp

Vacuum tap

Gas tap

OPTIONAL EQUIPMENT

Taps additional gaseous fluids.

Acces for different electrical equipment put on the car.

A second line of fuel gas

Spray formalin and related accessories to cycle semi-automatic decontamination.

Passive extraction kit for external expulsion.

Kit motorized bike extraction fan remote for external expulsion.

Tubular table support





IRUMABIO-

The last generation microbiological safety cabinets CRUMA BIO CLASS II Type A2 with digital functions.

The cabinet design ensures complete integrity by close control of the airflow. Incoming air is filtered through the first HEPA filter and down through the work area in a laminar flow pattern. The descending air creates a protective barrier preventing any outside air entering the cabinet. The airflow is then re-directed from the base of the cabinet into a plenum where 30% is expelled through a second HEPA filter and 70%, along with 30% new make up air, is re-circulated back into the cabinet through the first HEPA filter.



APPLICATIONS

- $\sqrt{}$ Manipulation of microorganisms, bacteria, fungi, viruses and parasites Risk categories 1, 2 and 3.
- \checkmark Isolation and sample culture
- √ Quantification methods
- $\sqrt{
 m Microscopy}$ techniques and sample preparation Identification and classification of microorganisms
- √ Genetic Manipulation

MAIN SPECIFICATIONS

- \checkmark Microprocessor controlled motor blower, with volumetric sensor for exhausted air flow monitoring
- √ State of the art Microprocessor control system offering:
 - Large screen monitor
 - Automatic control of preset airflow volumes
 - Sliding sash window with smart control
 - Permanent monitoring of HEPA filters life span
 - Alarms. Multilevel alarms, with redundancy functions Permanent display of working conditions
 - Highest air flow stability both in case of transitional disturbances or to progressive filter clogging
 - Semi-automatic fumigation cycle (EN12297 tested and certified)

More information www.cruma.es

- Continuous monitoring of front barrier air flow for the highest operator safety
- Low barrier alarm
- Power failure alarm
- √ Volt-free contact for remote monitoring of exhaust fan
- Automatic reset of initial conditions in case of powerfailure $\sqrt{}$
- \checkmark C-shaped support stand for the easiest one man installation procedure

CRUMABIO-E

CHARACTERISTICS

- √ Fully EN12469 certified by TÜV Hamburg
- \checkmark State of the art microprocessor control system
- Large digital display, high resolution $\sqrt{}$
- √ Air and aerosol-tight sliding sash, electrically operated by finger touch
- $\sqrt{}$ Alarms for low airflow and wrong front window position
- $\sqrt{}$ Sloped front and back wall for the most comfortable access
- $\sqrt{}$ Front access for filter maintenance and service
- C-shaped support stand for the easiest one man $\sqrt{}$ installation procedure
- √ Easy retrofit option kits

TECHNICAL SPECIFICATIONS					
Number of filters HEPA-H14		2			
Number of IP44 fans		1			
Average volume of treated air		±400 m³/h			
Average face velocity		0.50 m/s			
Total electrical power consumption		360 W			
Voltage-Frequency		110-230 V / 50-60 Hz			
Fluorescent Lamp / Light intensity		2x30W / 1200 Lux			
Noise level		50 dB			
Packaging:	Volume	2,39 m ³			
Palletized cardboard box	Weight	280 Kg			

2 year *warranty*

Because we are convinced of the quality of our products.





We recognise our responsability and dependence towards a healthy environment and, therefore, we destinate more than 7% of our annual budget in innovating and developing new products for the lab operator

our 72 engagement

Do you need help or technical assistance? Contact your distributor or call us if you have any questions or need technical support, spare parts,

maintenance service... +34 93 370 61 62

cruma.es export@cruma.es

(R) our D'engagement





SIZES (mm)					
External Dimensions		Internal Dimensions			
Widht 1380	Depth 840	Height 1450	Widht Depth Hei 1230 600 70		

MAIN STRUCTURE	
External metal parts	1.2 mm galvanized steel, coated with antiacid polymer powder resin thermo-hardened at 200 °C
Internal metal parts	Stainless Steel internal surfaces with 2B finishing (including spillage tray). Solid work surface in 3 sections and special designed front grill at the air barrier that cannot be obstructed by the operator's arms
Front door	Electrically operated from control panel sliding multilayer safety glass window

ACCORDING TO STANDARDS			
Cabinet	CEN-12469		
Filters	EN-779: 1996 (HEPA & ULPA Filters) EN-1822:1998 (HEPA & ULPA Filters)		
Quality	UNE EN ISO 9001:2015		



KEYBOARD MEMBRANE

Membrane keyboard commands for opening and closing the front window, activation of service outlets, electrical socket, gas solenoid valve, lighting and predisposed key to activate the connector kit germicidal lamp under the condition of complete closure of front of the machine, given the presence of the relevant position sensors.

All available accessories can be mounted without any structural modification of the cabin.

Provision for possible expulsion for channel adapters or additional filters.







MECHANICAL AND FUNCTIONAL SPECIFICATIONS

Sloped front design for the highest operational comfort. Sloped back side of the working chamber for the best down flow disbution (cabinet carcass EN12298 tested and certified for air tightness)

Utilities inlets from the top of the cabinet

Stainless Steel internal surfaces with 2B finishing (including spillage tray). Solid work surface (3 sections) and special designed front grill

Electrically operated sliding multilayer safety glass window (max opening at 120°)

Comfortable 200mm front opening

Easy to install retrofit options through lateral sides.

Exposed exhaust HEPA filter for easy visual integrity check

H14 class High Efficiency Particulate Airfilters with 99.995% efficiency on.3micron particles (most penetrating particle diameter) (EN 1822-1 and EN 13091:1999 tested and certified)

Both exhaust and Main Filters are equipped with a micromesh membran elocated downstream which acts as airspeed equalizer expansion plenum, as well as a clear indicator of filter damages

Filter change and maintenance from the front of the cabinet

Exhaust transitions easily installable

Key operated. The key can be removed when the unit is in SAFE mode, in order to avoid unwanted operation. In case of power failure, the cabinet is reset to original working conditions

Self calibration cycle performed when cabinet is switched on

High speed rinse and set up cycle performed, before reaching the SAFE operating mode.

Visual display of SAFE conditions. Pre-warning before actual alarm conditions are reached (visual and acoustic alarms)

Soft touch control with keys for standard service utilities. Interconnected UV and fluorescent lights

Exhaust and recirculating flow rates ensure 25 air changes/min in the working area (30% 70% split)

Front barrier air speed ≥0.5mt/sec

Aperture protection Factor (Apf)≥1.5 x 10 exp5

Cleanability Index CC grade. (EN 12296 tested and certified)

SERIAL EQUIPMENT

2 Electrical sockets

UV lamp

Vacuum tap

.

Gas tap

OPTIONAL EQUIPMENT

Taps additional gaseous fluids.

Acces for different electrical equipment put on the car.

A second line of fuel gas

Spray formalin and related accessories to cycle semi-automatic decontamination.

Passive extraction kit for external expulsion.

Kit motorized bike extraction fan remote for external expulsion.

Tubular table support





CRUMAVIRUS-E

Class III MICROBIOLOGICAL SAFETY CABINET

The recent emergence of Ebola virus in West Africa risks to be a growing worldwide threat!

Despite years of research on Ebola virus it is still not possible to deliver vaccines or treatments to the at-risk population or medical aid teams. Therefore there is currently no prophylaxis or treatment for Ebola virus infection.

This is the reason why the World Health Organization classifies Ebola virus as a pathogen of Risk Group Level 4. A Risk Group Level 4 pathogen must be manipulated, according to WHO, in a Containment Level 4 environment, wearing specific protective clothing and working with a Microbiological Safety Cabinet of Class III.

MAIN SPECIFICATIONS

- \checkmark Controls comfortably located at eye level
- \checkmark Fan speed and aeraulic controlled by Microprocessor
- \checkmark Three operating modes: normal, stand-by, calibration
- \checkmark High speed rinse at start up
- √ Self calibration and internal Watch-dog cycle before "SAFE" condition is reached
- √ Visual display of "SAFE" conditions and "UNSAFE" conditions (LED and bar graph)
- \checkmark Elapsed time meter
- \checkmark Microprocessor control with following specifications:
 - Multilevel alarms, with redundancy functions.
 - Permanent display of working conditions.
 - High air flow stability both in case of transitional disturbances or to progressive filter clogging
 - Power failure alarm
- \checkmark Volt-free contact for remote monitoring of exhaust fan.
- \checkmark Automatic reset of initial conditions in case of power failure
- \checkmark C-shaped support stand for easy one man installation procedure

More information www.cruma.es

- \checkmark Anti blow back valve (optional) for ducted configuration
- √ Magnehelic Gauge for internal chamber pressure constant monitoring
- \checkmark One (1) Electrical Socket as standard option
- \checkmark UV-Light installed on top (standard option)

TECHNICAL FEATURES LT20000 Order number N° of glove ports 2 Exhaust air flow rate (m3/h) >180 m3/h Internal Differential pressure (Pa) < -220 210 Weight (kg) Power Supply 220/240V 50Hz Power (W) 500 W Noise level < 58dB(A) Lighting lux >1000

Class III

SIZES (mm)						
External Dimensions			Internal Dimensions			
Widht 2105	Depth 822	Height 1300	Widht 1200	Depth 660	Height 700	

Do you need help or

technical assistance?

Contact your distributor or call us if you have any questions or need technical support, spare parts,

maintenance service...

CE (R) our D'engagement

+34 93 370 61 62

2 year *warranty*

Because we are convinced of the quality of our products.





cruma.es export@cruma.es

We recognise our responsability

and dependence towards a healthy environment and,

therefore, we destinate more than 7% of our annual budget in innovating and developing new products for the lab operator

our 70 engagemen

CRUMAURUS-2

TECHNICAL SPECIFICATIONS

- √ Manufactured in accordance with EN12469:2000 standard
- \checkmark State of the art microprocessor control system
- \checkmark Main switch with removable key
- \checkmark Soft touch keyboard
- ✓ Bar graph for exhaust air flow conditions; permanent display
- \checkmark Alarms for low air flow
- \checkmark Sloped front for the most comfortable access
- \checkmark Front and side access for filter maintenance and service
- \checkmark C-shaped support stand for easy one man installation procedure
- √ Transfer hatch with interlocked doors (right or left positioned)
- \checkmark Class III cabinet with exclusive three filter design and Class 100 inner chamber.





MECHANICAL AND FUNCTIONAL SPECIFICATIONS

- \checkmark Sloped front design for the highest operational comfort.
- \checkmark Stainless Steel internal surfaces with 2B finishing
- √ Liquid retaining work surface (Stainless Steel 2B finishing)
- ✓ Total visibility air and aerosol-tight front window equipped with robust gloves (Class III) for the safest operation when working with Risk Group 4 pathogens.
- ✓ Class III: Exclusive four filter design for the highest safety of the environment and the operator (Risk Group 4 pathogens): one (1) prefilter, one (1) HEPA H14 In-Let, two (2) HEPA H14 Exhaust Filters.
- ✓ H14 class High Efficiency Particulate Air filters with 99.995% efficiency on MPPS (most penetrating particle size) (EN1822-1 and EN 13091:1999 tested and certified)
- \checkmark Filter change and maintenance from the front of the cabinet.
- \checkmark Exhaust transitions easily installable.
- \checkmark Anti-blow-back valve available as add-on option
- ✓ Key operated. The key can be removed when the unit is in SAFE mode, in order to avoid unwanted operation. In case of power failure, the cabinet is re-set to original working conditions.
- \checkmark Self calibration cycle performed when cabinet is switched on.
- ✓ High speed rinse and set up cycle performed, before reaching the SAFE operating mode.
- √ Visual display of SAFE conditions. Pre-warning before actual alarm conditions are reached (visual and acoustic alarms)
- ✓ Soft touch control with keys for standard service utilities. Interconnected UV and fluorescent lights.
- √ 100% air exhaust.
- \checkmark Light intensity on work surface > 1000 lux.
- √ Noise level ≤58 dB(A)
- √ Work surface displacement (vibration) <0.005mm RMS between 20Hz and 20,000Hz (ISO 5349 tested and certified)
- √ 230V, 50Hz (230V,60Hz also available)
- \checkmark Max power (for each power point): 3Amps.
- \checkmark Microprocessor equipped with analogical watch-dog.
- \checkmark Leakage tested in agreement with EN 12469 and IS010648.2
- $\checkmark\,$ CE certification according to Machinery Directive 2006/42/CEE





CLASS III MICROBIOLOGICAL SAFETY CABINET





We recognise our responsability and dependence towards a healthy environment and, therefore, we destinate more than 7% of our annual budget in innovating and developing new products for the lab operator

our 7% engagement



cruma.es export@cruma.es





